

# EIF4A3 antibody

Catalog Number: 110258

## Product name

EIF4A3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

EIF4A3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, human lung tissue, mouse liver tissue, Raji cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 47kd

## Preparation

This antibody was obtained by immunization of EIF4A3 recombinant protein (Accession Number: NM\_014740). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110258(EIF4A3 antibody) at dilution of 1:1000



IP Result of anti-EIF4A3 (IP:Catalog No:110258, 3ug; Detection:Catalog No:110258 1:300) with HEK-293 cells lysate 2800ug.