

# EIF4H antibody



Catalog Number: 110265

## Product name

EIF4H antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

EIF4H Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 27kd

## Preparation

This antibody was obtained by immunization of EIF4H recombinant protein (Accession Number: NM\_031992). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110265(EIF4H antibody) at dilution of 1:800