

ELF1 antibody

Catalog Number: 110274

Product name

ELF1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

ELF1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human tonsillitis tissue. Positive WB detected in PC-3 cells, A431 cells, HeLa cells, HL-60 cells, Jurkat cells, mouse thymus tissue. Observed molecular weight by Western-blot: 97 kDa

Preparation

This antibody was obtained by immunization of ELF1 recombinant protein (Accession Number: NM_001370332). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110274(ELF1 antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human tonsillitis slide using Catalog No:110274(ELF1 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:110274(ELF1 Antibody) at dilution of 1:50