

# ELK1 antibody



Catalog Number: 110278

## Product name

ELK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ELK1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Jurkat cells, HepG2 cells, mouse pancreas tissue. Observed molecular weight by Western-blot: 62 kDa

## Preparation

This antibody was obtained by immunization of ELK1 recombinant protein (Accession Number: XM\_017029339). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110278(ELK1 antibody) at dilution of 1:800



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:110278(ELK1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)