

ELK4 antibody



Catalog Number: 110279

Product name

Polyclonal

ELK4 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

ELK4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, MCF7 cells, mouse liver tissue. Observed molecular weight by Western-blot: 47kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ELK4 recombinant protein (Accession Number: NM_021795). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110279(ELK4 antibody) at dilution of 1:800