

EBP2 antibody



Catalog Number: 110287

Product name

EBP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EBP2 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, A549 cells, HEK-293 cells, HeLa cells, NIH/3T3 cells, Raji cells. Observed molecular weight by Western-blot: 35-38 kDa

Preparation

This antibody was obtained by immunization of EBP2 recombinant protein (Accession Number: NM_006824). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110287(EBP2 Antibody) at dilution of 1:2000