

# ECE2 antibody

Catalog Number: 110291

## Product name

ECE2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ECE2 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse lung tissue, mouse ovary tissue. Observed molecular weight by Western-blot: 100 kDa

## Preparation

This antibody was obtained by immunization of ECE2 recombinant protein (Accession Number: BC005835). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110291(ECE2 antibody) at dilution of 1:500



IP Result of anti-ECE2 (IP:Catalog No:110291, 4ug; Detection:Catalog No:110291 1:500) with mouse liver tissue lysate 4000ug.