

EDRF1 antibody



Catalog Number: 110307

Product name

EDRF1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EDRF1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human placenta tissue. Positive WB detected in PC-3 cells, HEK-293 cells, HeLa cells, human brain tissue, MCF7 cells, mouse brain tissue. Observed molecular weight by Western-blot: 140-160kd

Preparation

This antibody was obtained by immunization of EDRF1 recombinant protein (Accession Number: NM_001202438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC

Dilutions


Recommended Dilution:

Validations




PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog

No:110307(C10orf137 antibody) at dilution of



1:1000

Immunohistochemical of paraffin-embedded human testis using Catalog No:110307(C10orf137 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:110307(C10orf137 antibody) at dilution of 1:50 (under 40x lens)