

ENPEP antibody

Catalog Number: 110322

Product name

ENPEP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ENPEP Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive FC detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of ENPEP recombinant protein (Accession Number: NM_001977). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, FC

Dilutions


Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using ENPEP antibody Catalog No:110322 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



1X10⁶ HepG2 cells were stained with 0.2ug ENPEP antibody (Catalog No:110322, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.