

ENPP5 antibody



Catalog Number: 110325

Product name

ENPP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ENPP5 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of ENPP5 recombinant protein (Accession Number: XM_011514786). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110325(ENPP5 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110325(ENPP5 antibody) at dilution of 1:100 (under 10x lens)