

# ENTPD4 antibody



Catalog Number: 110333

## Product name

ENTPD4 antibody

## Immunogen

[Human ENTPD4 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ENTPD4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human placenta tissue. Observed molecular weight by Western-blot: 70kd

## Preparation

This antibody was obtained by immunization of ENTPD4 recombinant protein (Accession Number: BC034477). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110333(ENTPD4 antibody) at dilution of 1:400