

EPHA6 antibody



Catalog Number: 110350

Product name

EPHA6 antibody

Specificity

Human; other species not tested.

Antibody description

EPHA6 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

Preparation

This antibody was obtained by immunization of EPHA6 recombinant protein (Accession Number: BC156733). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110350(EPHA6 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:110350(EPHA6 Antibody) at dilution of 1:200 (under 40x lens)