

# ERCC6 antibody

Catalog Number: 110352

## Product name

ERCC6 antibody

## Specificity

Human; other species not tested.

## Antibody description

ERCC6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 170 kDa

## Preparation

This antibody was obtained by immunization of ERCC6 recombinant protein (Accession Number: NM\_000124). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110352(ERCC6 Antibody) at dilution of 1:1000



IP Result of anti-ERCC6 (IP:Catalog No:110352, 4ug; Detection:Catalog No:110352 1:300) with HeLa cells lysate 3200ug.