

# ERCC6L antibody



Catalog Number: 110354

## Product name

ERCC6L antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

ERCC6L Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse lung tissue. Observed molecular weight by Western-blot: 180kd

## Preparation

This antibody was obtained by immunization of ERCC6L recombinant protein (Accession Number: BC008808). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110354(ERCC6L antibody) at dilution of 1:800