

ESCO2 antibody



Catalog Number: 110377

Product name

ESCO2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ESCO2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse testis tissue. Observed molecular weight by Western-blot: 65-70 kDa

Preparation

This antibody was obtained by immunization of ESCO2 recombinant protein (Accession Number: XM_011544421). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110377(ESCO2 antibody) at dilution of 1:1000