

# ESCO2 antibody



Catalog Number: 110377

## Product name

ESCO2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ESCO2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse testis tissue. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of ESCO2 recombinant protein (Accession Number: XM\_011544421). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110377(ESCO2 antibody) at dilution of 1:1000