

# Ephrin A4-Specific antibody



Catalog Number: 110385

## Product name

Ephrin A4-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Ephrin A4-Specific Rabbit Polyclonal antibody. Positive WB detected in SGC-7901 cells, human spleen tissue, mouse ovary tissue, rat lymph tissue. Observed molecular weight by Western-blot: 80kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_182690). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:110385(EFNA4-Specific antibody) at dilution of 1:500