

EPN3 antibody



Catalog Number: 110389

Product name

Polyclonal

EPN3 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

EPN3 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Dilutions

Preparation

Recommended Dilution:

This antibody was obtained by immunization of EPN3 recombinant protein (Accession Number: BC001038). Purification method: Antigen affinity purified.

WB: 1:500-1:5000

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Validations

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110389(EPN3 Antibody) at dilution of 1:1000