

ERBB4 antibody

Catalog Number: 110404

Product name

ERBB4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ERBB4 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 180 kDa

Preparation

This antibody was obtained by immunization of ERBB4 recombinant protein (Accession Number: NM_005235). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:500-1:5000

Validations



IP Result of anti-ERBB4 (IP:Catalog No:110404, 5ug; Detection:Catalog No:110404 1:1000) with mouse brain tissue lysate 3600ug.