

EXTL3 antibody



Catalog Number: 110419

Product name

EXTL3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EXTL3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of EXTL3 recombinant protein (Accession Number: XM_047421518). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110419(EXTL3 antibody) at dilution of 1:100



(under 25x lens)

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110419(EXTL3 Antibody) at dilution of 1:600