

EYA2 antibody



Catalog Number: 110421

Product name

EYA2 antibody

ELISA, WB, IHC

Specificity

Human, Mouse, Rat; other species not tested.

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Antibody description

EYA2 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in HeLa cells, mouse eye tissue, Y79 cells. Observed molecular weight by Western-blot: 57kd

Validations



Preparation

This antibody was obtained by immunization of EYA2 recombinant protein (Accession Number: BC013882). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110421(EYA2 antibody) at dilution of 1:100

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:110421(EYA2 antibody) at dilution of 1:50 (under 10x lens)