

EYA3 antibody

Catalog Number: 110422

Product name

EYA3 antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

EYA3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells. Positive IP detected in L02 cells. Positive IF detected in HeLa cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 63-66 kDa

Preparation

This antibody was obtained by immunization of EYA3 recombinant protein (Accession Number: NM_172098). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



Immunofluorescent analysis of HeLa cells, using EYA3 antibody Catalog No:110422 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110422(EYA3 Antibody) at dilution of 1:600



EYA3 antibody


Catalog Number: 110422

IP Result of anti-EYA3 (IP:Catalog No:110422, 4ug; Detection:Catalog No:110422 1:1000) with L02 cells lysate 1800ug.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

No:110422(EYA3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:110422(EYA3 Antibody) at dilution of 1:50 (under 40x lens)