

F2 antibody

Catalog Number: 110427

Product name

F2 antibody

Immunogen

[Human F2 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

F2 Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of F2 recombinant protein (Accession Number: NM_000506). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:5000-1:50000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:110427(F2 Antibody) at dilution of 1:30000