

F7 antibody



Catalog Number: 110431

Product name

F7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

F7 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human plasma tissue, mouse liver tissue. Observed molecular weight by Western-blot: 49-55 kDa

Preparation

This antibody was obtained by immunization of F7 recombinant protein (Accession Number: NM_019616). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110431(F7 Antibody) at dilution of 1:300