

# FAAH antibody

Catalog Number: 110435

## Product name

FAAH antibody

## Immunogen

[Human FAAH Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FAAH Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, A431 cells, rat testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 60-70 kDa

## Preparation

This antibody was obtained by immunization of FAAH recombinant protein (Accession Number: NM\_001441). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110435(FAAH Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110435(FAAH Antibody) at dilution of 1:50 (under 40x lens)