

FABP5 antibody

Catalog Number: 110441

Product name

FABP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FABP5 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer. Positive WB detected in human brain tissue, A375 cells, mouse skin tissue. Observed molecular weight by Western-blot: 15kd

Preparation

This antibody was obtained by immunization of FABP5 recombinant protein (Accession Number: NM_001444). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110441(FABP5 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:110441(FABP5 antibody) at dilution of 1:50 (under 10x lens)