

# Factor XII antibody



Catalog Number: 110445

## Product name

Factor XII antibody

Recommended Dilution:

WB: 1:200-1:2000

IHC: 1:20-1:200


## Immunogen

[Human F12 Recombinant protein \(His tag\)](#)

## Validations

## Specificity

Human, Mouse, Rat; other species not tested.




## Antibody description

Factor XII Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human colon cancer tissue. Positive WB detected in Jurkat cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 28 kDa

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110445(F12 antibody) at dilution of 1:500

## Preparation

This antibody was obtained by immunization of Factor XII recombinant protein (Accession Number: BC012390). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded human liver using Catalog No:110445(F12 antibody) at dilution of 1:50 (under 10x lens)

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality


Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC



## Dilutions

Immunohistochemical of paraffin-embedded

# Factor XII antibody



Catalog Number: 110445

---

human liver using Catalog No:110445(F12

antibody) at dilution of 1:50 (under 40x lens)