

FADS antibody

Catalog Number: 110448

Product name

FADS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FADS Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HepG2 cells, L02 cells, rat liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 50 kDa, 65-70 kDa

Preparation

This antibody was obtained by immunization of FADS recombinant protein (Accession Number: NM_001184891). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110448(FLAD1 Antibody) at dilution of 1:600



Immunofluorescent analysis of HepG2 cells using Catalog No:110448(FLAD1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)