

FAF1 antibody

Catalog Number: 110451

Product name

FAF1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAF1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, K-562 cells, mouse brain tissue, mouse testis tissue. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 66-74kd

Preparation

This antibody was obtained by immunization of FAF1 recombinant protein (Accession Number: NM_007051). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110451(FAF1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:110451(FAF1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog

FAF1 antibody



Catalog Number: 110451

No:110451(FAF1 antibody) at dilution of 1:50 (under 40x lens)