

FAM105A antibody

Catalog Number: 110457

Product name

FAM105A antibody

Specificity

Human; other species not tested.

Antibody description

FAM105A Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive WB detected in COLO 320 cells, human placenta tissue. Observed molecular weight by Western-blot: ~35 kDa

Preparation

This antibody was obtained by immunization of FAM105A recombinant protein (Accession Number: NM_019018). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110457(FAM105A Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:110457(FAM105A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog

FAM105A antibody



Catalog Number: 110457

No:110457(FAM105A Antibody) at dilution of 1:50 (under 40x lens)