

# FAM118B antibody



Catalog Number: 110463

## Product name

FAM118B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FAM118B Rabbit Polyclonal antibody. Positive IHC detected in human brain. Positive WB detected in HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 39kd

## Preparation

This antibody was obtained by immunization of FAM118B recombinant protein (Accession Number: NM\_024556). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110463(FAM118B antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:110463(FAM118B antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:110463(FAM118B antibody) at dilution of 1:100 (under 40x lens)