

FAM119A antibody



Catalog Number: 110464

Product name

FAM119A antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FAM119A Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive WB detected in L02 cells, HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of FAM119A recombinant protein (Accession Number: NM_001330134). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110464(FAM119A antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using FAM119A antibody Catalog No:110464 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).