

FAM127A antibody



Catalog Number: 110468

Product name

FAM127A antibody

Specificity

Human; other species not tested.

Antibody description

FAM127A Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of FAM127A recombinant protein (Accession Number: NM_001078171). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110468(FAM127A antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110468(FAM127A antibody) at dilution of 1:100 (under 25x lens)