

# ETS1 antibody



Catalog Number: 110486

## Product name

ETS1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ETS1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse thymus tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of ETS1 recombinant protein (Accession Number: BC017314). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110486(ETS1 Antibody) at dilution of 1:600