

# ETV6 antibody

Catalog Number: 110490

---

## Product name

ETV6 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ETV6 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HeLa cells. Observed molecular weight by Western-blot: 57kd

## Preparation

This antibody was obtained by immunization of ETV6 recombinant protein (Accession Number: NM\_001987). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:110490(ETV6 antibody) at dilution of 1:1000