

EVL antibody



Catalog Number: 110495

Product name

Polyclonal

EVL antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

EVL Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of EVL recombinant protein (Accession Number: NM_016337). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110495(EVL antibody) at dilution of 1:200