

Exo70 antibody

Catalog Number: 110500

Product name

Exo70 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Exo70 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, HEK-293 cells, mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 74 kDa

Preparation

This antibody was obtained by immunization of Exo70 recombinant protein (Accession Number: NM_001013839). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110500(EXoc7 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110500(EXoc7 antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-EXoc7 (IP:Catalog No:110500, 3ug; Detection:Catalog No:110500 1:600) with

Exo70 antibody



Catalog Number: 110500

mouse brain tissue lysate 2640ug.