

# EXOSC10 antibody



Catalog Number: 110508

## Product name

EXOSC10 antibody

## Specificity

Human; other species not tested.

## Antibody description

EXOSC10 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive IF detected in NIH/3T3 cells.

## Preparation

This antibody was obtained by immunization of EXOSC10 recombinant protein (Accession Number: NM\_001001998). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF


## Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110508(EXOSC10 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:110508(EXOSC10 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)