

FANCA antibody

Catalog Number: 110514

Product name

FANCA antibody

Specificity

Human; other species not tested.

Antibody description

FANCA Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 163 kDa

Preparation

This antibody was obtained by immunization of FANCA recombinant protein (Accession Number: NM_001018112). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110514(FANCA antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:110514(FANCA Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:110514(FANCA Antibody) at dilution of 1:100 (under 40x lens).