

FANCA-Specific antibody



Catalog Number: 110515

Product name

Polyclonal

FANCA-Specific antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

FANCA-Specific Rabbit Polyclonal antibody.
Positive WB detected in HEK-293 cells, HeLa cells.
Observed molecular weight by Western-blot:
163kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of
Peptide (Accession Number: NM_000135).
Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

HEK-293 cells were subjected to SDS PAGE
followed by western blot with Catalog
No:110515(FANCA-Specific antibody) at dilution
of 1:200