

FANCG antibody

Catalog Number: 110519

Product name

FANCG antibody

IP: 1:200-1:1000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FANCG Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in Jurkat cells, HeLa cells, K-562 cells, mouse testis tissue, mouse thymus tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 65-68 kDa

Preparation

This antibody was obtained by immunization of FANCG recombinant protein (Accession Number: NM_004629). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110519(FANCG antibody) at dilution of 1:100



IP Result of anti-FANCG (IP:Catalog No:110519, 4ug; Detection:Catalog No:110519 1:300) with mouse testis tissue lysate 4000ug.



Immunofluorescent analysis of HeLa cells using Catalog No:110519(FANCG Antibody) at dilution

FANCG antibody



Catalog Number: 110519

of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)