

FANCL antibody



Catalog Number: 110522

Product name

FANCL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FANCL Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 41-43kd

Preparation

This antibody was obtained by immunization of FANCL recombinant protein (Accession Number: NM_018062). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110522(FANCL antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110522(FANCL Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110522(FANCL Antibody) at dilution of 1:50 (under 40x lens)