

# FBF1 antibody

Catalog Number: 110536

## Product name

FBF1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FBF1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 69-70 kDa

## Preparation

This antibody was obtained by immunization of FBF1 recombinant protein (Accession Number: BC012332). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110536(FBF1 antibody) at dilution of 1:400



IP Result of anti-FBF1 (IP:Catalog No:110536, 4ug; Detection:Catalog No:110536 1:400) with mouse brain tissue lysate 5000ug.