

FBXO18 antibody



Catalog Number: 110548

Product name

Polyclonal

FBXO18 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

FBXO18 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 118-123kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of FBXO18 recombinant protein (Accession Number: BC032674). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110548(FBXO18 antibody) at dilution of 1:500