

FBXW5 antibody



Catalog Number: 110565

Product name

FBXW5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FBXW5 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, rat brain tissue. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of FBXW5 recombinant protein (Accession Number: NM_018998). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110565(FBXW5 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:110565(FBXW5 Antibody) at dilution of 1:50 (under 40x lens)