

FCN3 antibody

Catalog Number: 110570

Product name

FCN3 antibody

Specificity

Human; other species not tested.

Antibody description

FCN3 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, HepG2 cells, Y79 cells. Observed molecular weight by Western-blot: 35kd+50kd

Preparation

This antibody was obtained by immunization of FCN3 recombinant protein (Accession Number: NM_173452). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110570(FCN3 antibody) at dilution of 1:500