

# FAM46B antibody

Catalog Number: 110589

## Product name

FAM46B antibody

## Immunogen

[Human FAM46B Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

FAM46B Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of FAM46B recombinant protein (Accession Number: NM\_052943). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110589 (FAM46B Antibody) at dilution of 1:1000