

FAM65B antibody



Catalog Number: 110596

Product name

FAM65B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FAM65B Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 115kd,66kd

Preparation

This antibody was obtained by immunization of FAM65B recombinant protein (Accession Number: XM_017011527). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110596(FAM65B antibody) at dilution of 1:500