

# FAM72A antibody



Catalog Number: 110599

## Product name

Polyclonal

FAM72A antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

FAM72A Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 15-17 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of FAM72A recombinant protein (Accession Number: NM\_001385240). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:110599(FAM72A antibody) at dilution of 1:800