

FAM72A antibody



Catalog Number: 110599

Product name

FAM72A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAM72A Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 15-17 kDa

Preparation

This antibody was obtained by immunization of FAM72A recombinant protein (Accession Number: NM_001385240). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:110599(FAM72A antibody) at dilution of 1:800